

Improve patient safety
with Optiflow nasal high flow



Improve patient safety¹ during monitored anesthesia care

Nasal high flow (NHF) for airway management during monitored anesthesia care (MAC)

The aim of airway management during sedation is to maintain oxygenation. Drugs used for sedation and analgesia can cause respiratory depression and potential inadvertent periods of apnea.


Nasal high flow has been shown to:

- Reduce the risk of desaturation¹⁻⁴
- Reduce the incidence of airway related interventions^{1,2,5}
- Reduce the risk of procedural interruptions¹



Positive airway pressure increases end expiratory lung volume⁶⁻¹⁰

Nasal high flow significantly increased end-expiratory lung volume by 26% compared to low flow oxygen.⁶




Low flow vs Nasal high flow
Adapted from Corley et al., 2011

Patients at risk of desaturation


- Elderly
- Obstructive sleep apnea (diagnosed and undiagnosed)
- High Mallampati score
- Lung disease
- ASA Class ≥ 2
- High BMI

NHF facilitates CO₂ clearance¹¹⁻¹³


Cross-section of the upper airway immediately prior to inspiration compares CO₂ concentrations between natural breathing and Optiflow™ nasal high flow.**



Normal breathing



Optiflow NHF
(60 L/min)



** Computational Fluid Dynamic modeling. Fisher & Paykel Healthcare. Data on file.

Improve patient safety¹⁴ during airway management

Nasal high flow for airway management during induction of general anesthesia

The principle objective of airway management during anesthesia is to maintain adequate oxygenation.

Nasal high flow has been shown to:

- Provide effective preoxygenation¹⁴⁻¹⁶
- Provide short-term supplemental oxygenation^{14,17-19}
- Improve patient safety¹⁴



* For preoxygenation and short-term oxygen during intubation (for up to 10 minutes)

Patients at risk of desaturation or may be challenging to facemask ventilate

- High BMI
- Elderly
- Difficult airways
- Mask phobia
- Facial hair, e.g., beards
- Edentulous

ASA difficult airway recommendation²⁰

Optimizing oxygenation is recommended during:

- Preparation for difficult airway management
- Unanticipated and emergency difficult airway management
- Extubation of the difficult airway

The optimal solution for airway management



Heated humidification



Up to 70 L/min



Up to 100% O₂

Optiflow Filtered Nasal Interface with CO₂ Sampling

- Flexible sampling tube to detect exhaled CO₂ from either mouth or nose*

S M L

* For flow range applicable for this product, refer to the user instructions.



High-efficiency hydrophobic HEPA filter:



Minimum bacterial efficiency of at least 99.999%



Minimum viral efficiency of at least 99.996%



F&P Optiflow TRACE


Optiflow Filtered Nasal Interface with CO₂ sampling



Daily use Optiflow Oxygen Kit

Set up once at the start of the day and continue to use throughout


Can be used on up to 30 patients within a 24-hour period when combined with a filtered nasal interface



Anesthesia workstation mounting

Multiple mounting options available, based on the anesthesia workstation model for readily available access

In MAC, nasal high flow has been shown to provide effective oxygenation²¹⁻²³ and reduce the risk of desaturation¹⁻⁴ and airway related interventions^{1,2,5}




The use of nasal high flow reduced the incidence of hypoxia in patients undergoing elective gastroscopy with propofol sedation, compared to standard oxygen therapy.²

- ↓ 97% Decrease in jaw lifts
- ↓ 100% Decrease in desaturation events from 8.4% to 0%

Adapted from Lin et al., 2019.

Nay et al., 2021 | Moderate to high risk patients
The use of nasal high flow reduced the incidence of desaturation in patients undergoing GI endoscopy under deep sedation, compared to standard oxygen therapy.⁵

- ↓ 66% Decrease in airway maneuvers
- ↓ 72% Decrease in desaturation events




Adapted from Nay et al., 2021.

- ↓ 83% Decrease in risk of procedural interruptions
- ↓ 74% Decrease in risk of airway maneuvers
- ↓ 63% Decrease in risk of hypoxemia

Thirukenkatarajan et al., 2023.

Meta-analysis by Thirukenkatarajan et al. demonstrated the efficacy of Nasal High Flow (NHF) for reducing the risk of hypoxemia and requirement for airway maneuvers and procedural interruptions.

Forest plot comparing risk of hypoxemia between the NHF and control groups




*Mantel-Haenszel
NHF n=2024, Control n=2037.
Total events: HFNO = 115, Control = 371. Heterogeneity: Tau² = 0.48; Chi² = 58.05, df = 16 (P < 0.00001); I² = 72%.
Test for overall effect: Z = 4.62 (P < 0.00001). Test for subgroup differences: Chi² = 0.70, df = 3 (P = 0.87), I² = 0%



Favors NHF

Favors control

In general anesthesia, nasal high flow has been shown to improve patient safety¹⁴ by providing effective preoxygenation¹⁴⁻¹⁶ and extending safe apnea time^{14, 17-19}




Adapted from Hua et al., 2020




Hua et al., 2020 | Elderly patients

The use of nasal high flow significantly increased the proportion of elderly patients that reached an apnea time of 10 minutes, compared to standard oxygen therapy.¹⁸



Wong et al., 2019 | High BMI patients

The use of nasal high flow provided a longer safe apnea time and higher minimum SpO₂ in morbidly obese patients during anesthesia induction, compared to standard oxygen therapy.¹⁹



Adapted from Wong et al., 2019

Kuo et al., 2022.

Meta-analysis by Kuo et al. showed that Nasal High Flow (NHF) better enhanced arterial partial pressure of oxygen (PaO₂) and extended safe apnea time and is not inferior to facemask oxygenation in reducing the risk of desaturation during tracheal intubation.¹⁴

Forest plot of PaO₂ (mmHg) after pre-oxygenation between NHF and face mask oxygen (FMO) groups

Study or subgroup	Mean Difference IV, Random, 95% CI	Mean Difference IV, Random, 95% CI
Non-obese subjects		
Pillai 2016		13.00 [-26.60, 52.60]
Rajan 2018		12.00 [-69.25, 93.25]
Ng 2018		72.00 [21.44, 122.56]
Hua 2020		85.00 [34.54, 135.46]
Zhou 2021		112.00 [70.73, 153.27]
Osman 2021		129.00 [102.91, 155.09]
Lyons 2021		48.00 [-5.34, 101.34]
Subtotal (95% CI)		71.23 [33.07, 109.39]
Heterogeneity: Tau ² = 2029.26; Chi ² = 30.45, df = 6 (P < 0.0001); I ² = 80%		
Test for overall effect: Z = 3.66 (P = 0.0003)		
Obese subjects		
Heinrich 2014		46.00 [-24.68, 116.68]
Rosén 2021		37.00 [1.70, 72.30]
Subtotal (95% CI)		38.80 [7.22, 70.37]
Heterogeneity: Tau ² = 0.00; Chi ² = 0.05, df = 1 (P = 0.82); I ² = 0%		
Test for overall effect: Z = 2.41 (P = 0.02)		
TOTAL (95% CI)		64.86 [32.33, 97.40]
NHF n=192, FMO n=187		
Heterogeneity: Tau ² = 1838.40; Chi ² = 37.46, df = 8 (P < 0.00001); I ² = 79%		
Test for overall effect: Z = 3.91 (P < 0.0001)		
Test for subgroup differences: Chi ² = 1.65, df = 1 (P = 0.20); I ² = 39.3%		

Adapted from Kuo et al., 2022



Product Codes



Heaterbase:

MR810JHU
Heaterbase

1 each



Mounting:

AA501
DIN Rail Mounting Humidifier Pole
1 each

AA502
Dovetail Rail Mounting Humidifier Pole
1 each

AA503
Vertical Clamp Mounting Humidifier Pole
1 each



Oxygen Kits:

AA403
Optiflow Oxygen Kit - daily-use (for use with filtered nasal interface)
10/box

AA401SU
Optiflow Oxygen Kit - single use
10/box



Interfaces:

AA031J S / M / L
Optiflow Filtered Nasal Interface with CO₂ Sampling
10/box

AA030J S / M / L
Optiflow Nasal Interface with CO₂ Sampling
10/box

S M L



Transport:

AA020
Low Flow Oxygen Tubing
10/box

Support & Education

The F&P sales specialists provide a variety of education and support options for new and existing customers.



Scan to get in touch with our F&P sales representative

Optiflow Connect

Optiflow Connect allows new and existing users to access essential knowledge to apply Nasal High Flow (NHF) in their anesthesia practice.

Content includes:

- What NHF is and where it could be used.
- What patients could benefit from NHF.
- Clinical evidence for use during sedation, preoxygenation and during intubation.
- How to set up the equipment.

Clinician Website (secured access)

Clinician specific website is available by invite. It contains product and therapy information, along with additional clinical evidence and key opinion leader presentations.


Access to this site can be provided on request.

Education

Interactive educational sessions to provide opportunity for anesthesia and support staff to learn more about the use of NHF and discuss experiences.

Ongoing On-site Training

On-site training with F&P sales specialist, to provide education and product support onsite for new and existing users.



Scan the QR code to subscribe to our Quarterly Clinical Newsletter for recent research papers on nasal high flow in anesthesia.



References

- Thiruvenkatarajan V, Sekhar V, Wong DT, Currie J, Van Wijk R, Ludbrook GL. Effect of high-flow nasal oxygen on hypoxaemia during procedural sedation: a systematic review and meta-analysis. *Anaesthesia*. 2023;78(1):81-92.
- Lin Y, Zhang X, Li L, Wei M, Zhao B, Wang X, et al. High-flow nasal cannula oxygen therapy and hypoxia during gastroscopy with propofol sedation: a randomized multicenter clinical trial. *Gastrointest Endosc*. 2019;90(4):591-601.
- Hung KC, Chang YJ, Chen IW, Soong TC, Ho CN, Hsing CH, et al. Efficacy of high flow nasal oxygenation against hypoxemia in sedated patients receiving gastrointestinal endoscopic procedures: A systematic review and meta-analysis. *J Clin Anesth*. 2022;77:110651.
- Su CL, Chiang LL, Tam KW, Chen TT, Hu MC. High-flow nasal cannula for reducing hypoxic events in patients undergoing bronchoscopy: A systematic review and meta-analysis of randomized trials. *PLoS One*. 2021;16(12):e0260716.
- Nay MA, Fromont L, Eugene A, Marcueyz JL, Mfam WS, Baert O, et al. High-flow nasal oxygenation or standard oxygenation for gastrointestinal endoscopy with sedation in patients at risk of hypoxaemia: a multicentre randomised controlled trial (ODEPHI trial). *Br J Anaesth*. 2021;127(1):133-42.
- Corley A, Caruana LR, Barnett AG, Tronstad O, Fraser JF. Oxygen delivery through high-flow nasal cannulae increase end-expiratory lung volume and reduce respiratory rate in post-cardiac surgical patients. *Br J Anaesth*. 2011;107(6):998-1004.
- Parke RL, Bloch A, McGuinness SP. Effect of Very-High-Flow Nasal Therapy on Airway Pressure and End-Expiratory Lung Impedance in Healthy Volunteers. *Respir Care*. 2015;60(10):1397-403.
- Parke RL, Eccleston ML, McGuinness SP. The effects of flow on airway pressure during nasal high-flow oxygen therapy. *Respir Care*. 2011;56(8):1151-5.
- Yuan Z, Han X, Wang L, Xue P, Sun Y, Frerichs I, et al. Oxygen Therapy Delivery and Body Position Effects Measured With Electrical Impedance Tomography. *Respir Care*. 2020;65(3):281-7.
- Tatsuiishi W, Sato T, Kataoka G, Sato A, Asano R, Nakano K. High-Flow Nasal Cannula Therapy With Early Extubation for Subjects Undergoing Off-Pump Coronary Artery Bypass Graft Surgery. *Respir Care*. 2020;65(2):183-90.
- Möller W, Celik G, Feng S, Bartenstein P, Meyer G, Eickelberg O et al. Nasal high flow clears anatomical dead space in upper airway models. *J Appl Physiol*. 2015;118(12):1525-32.
- Pinkham M, Burgess R, Mündel T, Tatkov S. Nasal high flow reduces minute ventilation during sleep through a decrease of carbon dioxide rebreathing. *J Appl Physiol*. 2019;126(4):863-9.
- Onodera Y, Akimoto R, Suzuki H, Nakane M, Kawamae K. A high-flow nasal cannula system set at relatively low flow effectively washes out CO(2) from the anatomical dead space of a respiratory-system model. *Korean J Anesthesiol*. 2017;70(1):105-6.
- Kuo HC, Liu WC, Li CC, Cherng YG, Chen JT, Wu HL, et al. A comparison of high-flow nasal cannula and standard facemask as pre-oxygenation technique for general anesthesia: A PRISMA-compliant systemic review and meta-analysis. *Medicine*. 2022;101(10):e28903.
- Heinrich S, Horbach T, Stubner B, Prottengeier J, Irouschek A, Schmidt J. Benefits of heated and humidified high flow nasal oxygen for preoxygenation in morbidly obese patients undergoing bariatric surgery: a randomized controlled study. *J Obes Bariatrics*. 2014;1(1):7.
- Zhou S, Zhou Y, Cao X, Ni X, Du W, Xu Z, et al. The efficacy of high flow nasal oxygenation for maintaining maternal oxygenation during rapid sequence induction in pregnancy: A prospective randomised clinical trial. *Eur J Anaesthesiol*. 2021;38(10):1052-8.
- Osman YM, Abd El-Raof R. High flow nasal cannula oxygen preventing desaturation during induction of general anaesthesia in caesarean section: A randomized controlled trial. *Trends Anaesth Crit Care*. 2021;40:23-7.
- Hua Z, Liu Z, Li Y, Zhang H, Yang M, Zuo M. Transnasal humidified rapid insufflation ventilatory exchange vs. facemask oxygenation in elderly patients undergoing general anaesthesia: a randomized controlled trial. *Sci Rep*. 2020;10(1):5745.
- Wong DT, Dallaire A, Singh KP, Madhusudan P, Jackson T, Singh M, et al. High-Flow Nasal Oxygen Improves Safe Apnea Time in Morbidly Obese Patients Undergoing General Anesthesia: A Randomized Controlled Trial. *Anesth Analg*. 2019;129(4):1130-6.
- Apfelbaum JL, Hagberg CA, Connis RT, Abdelmalak BB, Agarkar M, Dutton RP, et al. 2022 American Society of Anesthesiologists Practice Guidelines for Management of the Difficult Airway. *Anesthesiology*. 2022;136(1):31-81.
- Kim SH, Bang S, Lee KY, Park SW, Park JY, Lee HS, et al. Comparison of high flow nasal oxygen and conventional nasal cannula during gastrointestinal endoscopic sedation in the prone position: a randomized trial. *Can J Anaesth*. 2021;68(4):460-6.
- Irfan M, Ahmed M, Breen D. Assessment of High Flow Nasal Cannula Oxygenation in Endobronchial Ultrasound Bronchoscopy: A Randomized Controlled Trial. *J Bronchology Interv Pulmonol*. 2021;28(2):130-7.
- Lee MJ, Cha B, Park JS, Kim JS, Cho SY, Han JH, et al. Impact of High-Flow Nasal Cannula Oxygenation on the Prevention of Hypoxia During Endoscopic Retrograde Cholangiopancreatography in Elderly Patients: A Randomized Clinical Trial. *Dig Dis Sci*. 2022;67(8):4154-60.
- Ferrando C, Puig J, Serralta F, Carrizo J, Pozo N, Arocas B, et al. High-flow nasal cannula oxygenation reduces postoperative hypoxemia in morbidly obese patients: a randomized controlled trial. *Minerva Anestesiol*. 2019;85(10):1062-70.
- Shih CC, Liang PC, Chuang YH, Huang YJ, Lin PJ, Wu CY. Effects of high-flow nasal oxygen during prolonged deep sedation on postprocedural atelectasis: A randomised controlled trial. *Eur J Anaesthesiol*. 2020;37(11):1025-31.
- Xu L, Wei N, Shan M, Wang Z, Ye T, Liang S, et al. [High-flow nasal cannula oxygen therapy can reduce occurrence of hypoxia in elderly patients during anesthesia recovery]. *Nan Fang Yi Ke Da Xue Xue Bao*. 2021;41(8):1265-9.

F&P, Optiflow, Optiflow Trace and THRIVE are trademarks of Fisher & Paykel Healthcare. For patent information, please see www.fphcare.com/ip. 614688 REV G © 2024 Fisher & Paykel Healthcare Limited.