

# THINK OUTSIDE THE STENT



**Medtronic**  
Further, Together



**amg**international  
A Q3 Medical Company

## **ARCHIMEDES™\***


biodegradable biliary and pancreatic stent

The ARCHIMEDES™\* biodegradable biliary and pancreatic stent is intended to be used to drain obstructed biliary or pancreatic ducts.<sup>1</sup> The patented double helical design of the stent allows bile or pancreatic fluid to flow on the outer extremity of the device while supporting the patency of the duct.

**A new way of thinking about biliary drainage —  
procedural intervention should no longer be the norm.**

# ARCHIMEDES™\*

biodegradable biliary and pancreatic stent



## Fast degrading profile\*

### 2 mm diameter (6 F)

Product code	Length (mm)
BPS20040F	40
BPS20060F	60
BPS20080F	80
BPS20100F	100
BPS20125F	125

### 2.6 mm diameter (~8 F)

Product code	Length (mm)
BPS26040F	40
BPS26060F	60
BPS26080F	80
BPS26100F	100
BPS26125F	125

### 3.4 mm diameter (~10 F)

Product code	Length (mm)
BPS34040F	40
BPS34060F	60
BPS34080F	80
BPS34100F	100
BPS34125F	125

## Slow degrading profile\*

### 2 mm diameter (6 F)

Product code	Length (mm)
BPS20040S	40
BPS20060S	60
BPS20080S	80
BPS20100S	100
BPS20125S	125

### 2.6 mm diameter (~8 F)

Product code	Length (mm)
BPS26040S	40
BPS26060S	60
BPS26080S	80
BPS26100S	100
BPS26125S	125

### 3.4 mm diameter (~10 F)

Product code	Length (mm)
BPS34040S	40
BPS34060S	60
BPS34080S	80
BPS34100S	100
BPS34125S	125

\* The suitable degradation profile of the stent to treat the obstructed biliary or pancreatic duct must be chosen by a clinical professional, always taking the underlying disease and the condition of the individual patient into account.

## Intended use/indications

This device is used to drain obstructed biliary or pancreatic ducts and is indicated for obstructed biliary or pancreatic ducts.

## Contraindications

- Those specific to ERCP and any procedure to be performed in conjunction with stent placement.
- Inability to pass the guidewire or stent through the obstructed area.
- All cardiovascular applications.
- Placement of the stent in biliary obstructions precluding any form of cholangiography.

## Reference

1. ARCHIMEDES instructions for use. (IFU ARCHIMEDES amg. 02 - 093A45\_180706)

©2020 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All rights reserved. Medtronic and the Medtronic logo are trademarks of Medtronic. ™\* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. 20-weu-archimedes-sell-sheet-bps-4032960

**Medtronic**  
Further, Together



**amg**international  
A Q3 Medical Company