




# Biliary Slim™ M Stent

6Fr delivery system with tapered tip



- Radiopaque marker: 3 (three) at both ends & 2 (two) in the middle
- Hook and Cross for conformability
- Delivery system : 6Fr / 180cm (0.025 inch guidewire)
- Stent size: Diameter 6,8mm / Length 4~12cm

TaeWoong Medical USA  
4801 Wilshire Blvd. Suite 280, Los Angeles, CA 90010  
E: [info@taewoongusa.com](mailto:info@taewoongusa.com) T: 213-944-0503  
[www.taewoongusa.com](http://www.taewoongusa.com)

©2025 TaeWoong Medical USA. All Rights Reserved.



## Features




### Versatility Excelling in Both SIS and SBS techniques

For MHBO cases requiring long-term drainage, the combination of low shortening rates, minimal tissue damage, and improved stability makes Slim M a strong choice for clinicians.



### Shinkansen Tip Advantage

Minimal guidewire-tip gap, enhanced passage through strictures and sharp angles, soft design for difficult cases.



### Minimized Tissue Damage with Soft Braided Structure

Slim M's braided structure minimizes unnecessary tissue irritation. The M-Type design reduces stent-induced stress, improving procedural outcomes and minimizing adverse events.

## Ordering Information

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (fr)	Usable Length (cm)*
BN0604-6	6	4	6	180
BN0605-6		5		
BN0606-6		6		
BN0607-6		7		
BN0608-6		8		
BN0609-6		9		
BN0610-6		10		
BN0612-6		12		
BN0804-6	8	4	6	180
BN0805-6		5		
BN0806-6		6		
BN0807-6		7		
BN0808-6		8		
BN0809-6		9		
BN0810-6		10		
BN0812-6		12		

### Indications for Use:

The Niti-S Biliary Slim M Stent is indicated for the palliation of malignant strictures in the biliary tree.