


FiAPC[®] plus

Precision starts with control


- Enhanced tissue effect
- Improved usability
- Electrosurgical expertise



Precise beam control with homogeneous tissue effect




Improved thermals



Increased flexibility



Improved markings




Improved suction



Ergonomic plug design

Technical Data

Tip configurations



AXIAL



20132-321 FiAPC® plus A; Ø 2.3mm; L 2.2m;
flexible probe w/ filter

20132-323 FiAPC® plus A; Ø 2.3mm; L 3.0m;
flexible probe w/ filter

CIRCUMFERENTIAL



20132-325 FiAPC® plus C; Ø 2.3mm; L 2.2m;
flexible probe w/ filter

SIDE FIRE



20132-324 FiAPC® plus SC; Ø 2.3mm; L 2.2m;
flexible probe w/ filter